

PORTABLE X-RAY UNIT

MIKASA

HF100H



HIGH FREQUENCY
INVERTER SYSTEM



MIKASA
Established in 1923

*HIGH POWER!
&
LIGHT WEIGHT!*

MIKASA PORTABLE X-RAY SERIES




MIKASA HF8015H



MIKASA HF100H

	HF8015H	HF100H
Output	15mA@50-80kV	40mA@40-44kV 30mA@46-60kV 25mA@62-80kV 20mA@82-100kV
Timer	0.01-2.00 sec.	0.02-8.00 sec.
AC input	230V(50/60Hz) specify voltage desired	230V(50/60Hz) specify voltage desired
Focal spot	0.8mm	1.2mm
Total filtration	2.7mm Al equivalent	2.7mm Al equivalent
Memory storage	5 presets	5 presets
Size	18.2cmx15.7cmx29.9cm 7.2"Wx6.2"Hx11.8"L	28.6cmx23.0cmx40.0cm 11.3"Wx9.1"Hx15.8"L
Weight	15.1 lbs (6.85 kgs)	37.0 lbs (16.8 kgs)



Control panel

5 - 10 - 5, Koishikawa, Bunkyo - ku, Tokyo 112 - 0002 JAPAN

Phone: 81 - 3 - 3811 - 6958 Fax: 81 - 3 - 3814 - 7844

E - mail: info@mikasax - ray.com <http://www.mikasax - ray.com>

MIKASA X-RAY CO., 