

PORTABLE X-RAY UNIT

MIKASA

HF400V / RHF-1B

HIGH FREQUENCY INVERTER SYSTEM

FOR VETERINARY USE

High output model



MIKASA High Frequency Portable X-ray unit, offer the highest power to weight ratio. The HF400V and RHF-1B are high output models for the equine and the large animal practitioner. MIKASA equipment is capable of all routine radiographic views you would expect from a portable unit, and compatible with all film and digital imaging systems.



CERTIFICATE ISO 9001 ISO 13485



MIKASA
Established in 1923

**HIGH POWER
&
LIGHT WEIGHT**



MIKASA PORTABLE X-RAY SERIES *High output models*




MIKASA HF400V



MIKASA RHF-1B

- ◆ Max. 100mA / Max. 120kV high output model...HF400V
- ◆ Powerful and light weight (11.2kgs) for the general Veterinary...RHF-1B
- ◆ Integrated, adjustable dual laser pointers
- ◆ 5 memory stations for storing x-ray exposure techniques
- ◆ Interface connector for DR system

	HF400V	RHF-1B
Max. output power	6.0 kW	2.5 kW
kV & mA range	(0.01-0.1sec.) 100mA@40-50kV 80mA@52-60kV 70mA@62-80kV 60mA@82-100kV 40mA@102-120kV (0.11-2.0sec.) mA 30% down	40mA@30-60kV 35mA@62-70kV 30mA@72-80kV 25mA@82-100kV
Max. mAs	143 mAs	60 mAs
AC input	230V(50/60Hz)	230V(50/60Hz)
Focal spot	1.2mm	1.2mm
Total filtration	2.7mm Al equivalent	2.7mm Al equivalent
Size	W340mm x H250mm x L430mm	W240mm x H196mm x L365mm
Weight	21.8 kgs	11.2 kgs



*Specifications are subject to change without notice.

5-10-5, Koishikawa, Bunkyo-ku, Tokyo 112-0002 JAPAN

Phone: 81-3-3811-6958 Fax: 81-3-3814-7844

E-mail: info@mikasax-ray.com <http://www.mikasax-ray.com>

MIKASA X-RAY CO.,LTD